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Operating Manual: BoneTrust® Implant Systems

Instructions for Use BoneTrust® implant systems

SCOPE OF THESE INSTRUCTIONS FOR USE

BoneTrust® implant systems and prosthetics in accordance with the Medical Device Directive 93/42/EEC.

PRODUCT INFORMATION

General description

The BoneTrust® implant systems manufactured by Medical Instinct® Deutschland GmbH comprise implants, transmucosal abutments and other prosthetic components, as well as surgical and prosthetic instruments. The components of the BoneTrust® implant systems are parts of an integrated system and should only be used in conjunction with other original components and original instruments.

Intended use and indications

Implants

BoneTrust® plus/sinus/cone+/cone+ sinus/cone PWR / hex PWR/one implants can be used to treat partial or total tooth loss. The implants support both single tooth restorations and provisionally detachable or fixed and removable bridgework and prostheses.

BoneTrust® mini/mini+ implants are interim implants exclusively for the provisional treatment of bridges or dental prostheses.

With BoneTrust® plus/sinus/cone+/cone+ sinus/cone PWR/hex PWR/ one implants, subgingival healing is recommended, whereby the treating physician may also opt for transgingival implantation by means of healing abutments (gingiva formers). BoneTrust® sinus implants are used for standardised implant treatment in the area of the maxillary sinus. The treating physician can opt for either manual or machine mounting.

The BoneTrust® implants with a reduced diameter (\varnothing 3.0 mm) are limited in terms of their indication and are intended for use in patients with a small alveolar ridge (at least 5.0 mm) with low stress. Owing to the lower mechanical strength compared to the implants with a diameter of 3.4 mm and 4mm, these implants are only suitable for the following indications:

- single tooth replacement for the lateral incisors in the maxilla and mandible, and of the central incisors in the mandible
- fixed dental prosthesis, in blocked prosthetic constructions and in conjunction with implants with a standard diameter of 4.0 mm
- removable dental prosthesis, with at least 4 implants in primary blocking
- as supplementary abutment

BoneTrust® mini/mini+ implants are interim implants exclusively for the provisional treatment of bridges and dental prostheses. The implants should be removed after at most 9 months.

Prosthetics

The system includes various different superstructures used to anchor the prosthetic implant, which are individually selected to match the user's indication and anatomical requirements. The prosthetics for the BoneTrust® implant systems manufactured by Medical Instinct® Deutschland GmbH offer a range for all prosthetic indications such as:

- single tooth and bridgework (cement-retained and occlusally screwed)
- hybrid constructions
- temporary restoration (single tooth and splinted)
- bar restoration
- occlusally screwed hybrid prostheses
- telescopic restorations
- prosthesis fixation

Customized threaded connections or interlocking geometries (hexagon and conical torx) ensure traction between interfaces and the individual combination devices.

Contraindications

Implants

There are no absolute contraindications for the use of BoneTrust® dental implants apart from those generally applicable to implant surgery.

The following situations may reduce the chances of successful implant treatment:

- General disorders and pathological maxillary diseases.
- Impaired healing capacity and malnutrition due to disease or (radiation) therapy and due to alcohol/drug abuse.
- Smoking also affects the outcome of the implantation.

Although the materials used are biocompatible materials, there is a low risk common to all types of implant dentistry that implants will not integrate into the bone. Implants or bone grafts may occasionally fail due to overloading of the bone or bacterial infections; in rare cases this may cause damage to sensory nerves or adjacent dental roots. Risks include poorly controlled diabetes with fluctuating blood sugar levels, poor oral hygiene and osteoporosis. There is an increased risk of implant failure during oral and/or intravenous administration of bisphosphonates. Absolute contraindications are infected extraction sockets, severe apical osteitis (bone inflammation) and bony defects.

Prosthetics

Use of incompatible instruments or combination devices.

Risks and side effects

The risks, complications and side effects may be intra-operative or post-operative.

Intra-operative risks, side effects and complications

Damage to major structures:

- inferior alveolar nerve (nervus alveolaris inferior);
- lingual nerve (nervus lingualis);
- adjacent teeth;
- soft tissue;
- maxillary sinus;
- hemorrhaging (vascular).

Post-operative risks, side effects and complications

Early-onset complications:

- secondary bleeding/hematoma
- edema;
- infection;
- wound healing impairment.

late-onset complications:

- periimplantitis;
- implant failure.

SAFETY

General

For dental use by specialist personnel only.

Read these Instructions for Use in full. Retain the Instructions for Use for reference. Failure to follow the Instructions for Use may lead to injuries or damage to the device.

Prior to the start of treatment each patient must undergo diligent medical and dental assessment. Prior to the implant treatment any pathological diseases of the jaw and remaining dentition must be fully healed, or the healing process must have begun. Bone quality/quantity, local infections or impairment of the initial healing process are factors that may affect integration of the implant. In order to determine the optimum implant length and to prevent damage to the nerve, the available space must be clarified in advance (x-ray). Application and selection of appropriate implants should be carried out by surgically trained medical specialists only.

Safety instructions

NOTE

NOTES are general precautions, which, if not followed, may lead to impairment or temporary inconvenience.

The word CAUTION is used to signal hazards which, if not followed, may lead to mild to moderate injury or impairment of the treatment.

The word WARNING is used to signal hazards which, if not followed, may lead to severe injury or severe impairment of the treatment.

CAUTION

WARNING

General safety instructions, liability and warranty

CAUTION

Use of third-party components may impair the function of the BoneTrust® implant systems and invalidates any warranty or liability for damages by Medical Instinct® Deutschland GmbH.

This applies particularly to application techniques other than those recommended. The processing and application of Medical Instinct® Deutschland GmbH devices are conducted outside of our control and are the sole responsibility of the user. No liability is accepted for damage so caused.

CAUTION

Take a medical history prior to each treatment and explain the contra-indications and potential adverse effects to the patient.

WARNING

When applied intraorally, ensure that measures are taken to prevent devices from being aspirated or falling out.

WARNING

Prior to application, check that the device is suitable for use in accordance with its intended purpose and indication. Check the color coding and lettering to prevent confusion over device models.

WARNING

The implants are single-use devices and are supplied sterile. Re-use or re-sterilization of the implants is not permitted. The implants are sterilized by gamma radiation. Re-using single-use devices increases the risk of infection and risk-free functional safety cannot be guaranteed.

WARNING

Inspect the devices and the inner packaging for damage prior to use. Devices that are damaged devices or arrive with damaged packaging must not be used. Sterility is not guaranteed if packaging is damaged or opened. Use of the devices after the expiration date is not permitted.

NOTE

The user is required to report all serious incidents occurring in connection with the product to the manufacturer and the responsible authority of the state in which the user is established.

USE

General preparations

WARNING

Inspect the device for damage prior to use. Damage may lead to injury and may impair the effectiveness of the treatment.

The prosthetic device is supplied non-sterile and must therefore be pre-sterilized. See the "Pre-processing" chapter.

Treatment plan

WARNING

Deficiencies in the patient's medical history, pre-operative diagnostics and treatment plan may lead to implant failure at an early stage.

Clinical and radiological examination of the patient, as well as a model analysis, are essential prerequisites for successful implant treatment. Conduct pre-operative diagnostics carefully and document the treatment plan.

Procedure for implant placement

WARNING

Primary implant stability is essential after the implant is inserted. Instruments and tools must be secured in order to prevent them being inadvertently aspirated by the patient.

Treatment must be hard and soft tissue friendly to ensure optimum conditions are met for integration of the BoneTrust® implant systems. The intervention should be conducted as atraumatically as possible and requires a high degree of precision and care. Parallel drilling with suitably sharp drills should, depending on the drill model, be carried out at 300-1,000 rpm (cf. Instructions for Use for BoneTrust® Implantology Drill Set), ensuring the drill tip is adequately cooled to avoid thermal trauma, which can prevent successful osseointegration of the implant. We recommend using physiologic, sterile saline solution as coolant.

NOTE

For the instructions for use, guides and documents, make sure that they have a current revision.

CAUTION

Please refer to the relevant Instructions for Use (Instructions for Use for BoneTrust® Implantology Drill Set) for further details on using compatible placement tools.

Placement of the implants, both machine and manual, should be performed at a max. torque of 35 Ncm (25 rpm).

CAUTION

The maximum torques to be used for BoneTrust® mini/+ implants are:
maximum insertion torque = 35 Ncm
maximum removal torque = 80 Ncm

Surgical procedure

WARNING

Any contamination of the implant must be avoided.

Implantation is a surgical procedure. It must be performed in accordance with the rules of surgery (explanation, sterility, aftercare). It is essential that treatment is hard and soft tissue friendly to ensure optimum osseointegration conditions are met.

1. Prepare the implant site in accordance with the selected surgical procedure and implant size and length. You should follow the recommended surgical guidelines (Surgical Guidelines for BoneTrust® mini and mini+ Interim Implants, 08/2020 Version and Surgical Guidelines for BoneTrust® Implant Systems in Device Catalog, 08/2020 Version pp. 16/17, Surgical Guidelines for BoneTrust PWR).
2. Remove the outer paper packaging, open the blister pack and carefully take out the glass container.
3. Remove the implant unit from the glass tube by pulling gently on the stopper.
4. Next, place the implant into the prepared implant side by rotating the stopper gently and then, once the implant is held in the drill hole, remove the stopper by pulling the mount or implant head (BoneTrust® and BoneTrust® mini).
5. Next, engage the corresponding placement tool fully onto the hex of the implant mount (implant head in the case of BoneTrust® mini/mini+ implants). Screw in the implant at low speed (max 25 rpm/35 Ncm). The ratchet may be used for final positioning to ensure exact alignment.
6. Labial/buccal alignment of one corner is necessary to allow optimal positioning of the superstructures later on.
7. Epicrestal or slightly subcrestal placement of the implant is recommended.
8. The maximum recommended torque for implant placement is 35 Ncm (25 rpm), regardless of implant diameter.

Wound closure

The implant must be locked in place before the gingiva is sutured post-placement. This may be accomplished using a screw plug or healing abutment.
Inspect the inside of the implant for any blood or tissue residues. Any residues visible inside the implant must be removed before the implant is locked in place. The retaining screw plug and the healing abutment should be screwed gently by hand only.

Next, complete a saliva-proof suture that is free of tension. The suture may be removed after 8 to 10 days with normally healing of the soft tissue. Prompt surgical revision with coverage is required if there is potential onset of wound dehiscence.

Healing phase

The mandibular healing phase normally lasts three months, the maxillary phase six months. It may be shorter or longer than this depending on bone quality and anatomy. Once the healing phase and gingival formation are complete, prosthetic restoration can start.

Attaching the superstructure

1. Place the selected abutment into the implant. Ensure that the abutment is positioned correctly. There should be no visible gap between the superstructure and the implant. If a visual check is not possible, check that it is seated correctly using a probe or if necessary by x-ray.
2. Connect the screwdriver to the torque ratchet and introduce the retaining screw into the designated hole on the abutment.
3. You can then attach the abutment to the implant with the retaining screw. The recommended torque for superstructure retaining screws is 20 Ncm.
4. Ensure that the screwdriver is located exactly in the screw head. Worn, third-party or screwdrivers that have been tampered with will damage the screw when it is tightened. A new, unused retaining screw should always be used for the final placement to help prevent it from loosening later on.
5. Alternatively, the retaining screw may also be placed using a handpiece. To do this, attach the screwdriver to the handpiece. The retaining screw is placed in the screw channel and tightened with the handpiece (max. 20 Ncm).

CAUTION

Check the precise torque is used when the handpiece is used to screw the mount in place. If your handpiece does not support torque measurement we recommend screwing manually with the torque ratchet.

Once the mount is screwed in place the screwdriver can be detached from the retaining screw by pulling gently.

Post-processing/final meeting

Inform the patient in detail on the oral hygiene measures required. Recommend that the patient should attend follow-up examinations regularly.

NOTE

Before the final prosthetic implantation the attachment screw should be tightened once more. Ideally, the attachment screw on the prosthetic implant should be given another final tightening after 3 months, wherever possible.

INSTRUCTIONS ON CLEANING AND STERILIZATION

Bone Trust® implant system implants are supplied sterile. They must not be re-processed. Medical Instinct® Deutschland GmbH accepts no liability for implants sterilized by the user.

CAUTION

The prosthetics and impression components are intended for single use only and like the placement tools must be prepared before first use since they are supplied non-sterile.

A validated reprocessing procedure must be used. Only disinfectants that have been tested, approved (e.g. VAH or FDA approval or CE marking) and are suitable for the instruments may be used for reprocessing the instruments. The manufacturer's instructions for use must be observed.

Thermal disinfection

When selecting the disinfectant ensure that it has proven efficacy (e.g. DGHM or FDA approved, or CE Mark to DIN EN ISO 15883). Thermal disinfection should be conducted prior to use on the patient using the corresponding small parts holder.

CAUTION

If any visible soiling is detected after machine processing, cleaning and disinfection should be repeated until soiling is no longer visible.

The cleaned components should be dried and sealed into sterilization pouches (in accordance with DIN EN 868-5).

Sterilization

Sterilization must be conducted using the pulsed flow vacuum process (Class B small steam sterilizer in accordance with DIN EN 13060). The minimum sterilization hold time (period of exposure to sterilization temperature) is 5 minutes at a minimum temperature of 134°C (273°F). Minimum drying time is 10 minutes. Processing of medical devices ends with the documented release for use.

CAUTION

The devices should be classed as "not ready for use" until processing is complete.

INFORMATION ON MRI SAFETY

Titanium is viewed as the gold standard in dental implantology and has been exhaustively studied in the relevant literature. When Magnetic Resonance Imaging is used during examinations and treatment, please note that the BoneTrust® implant systems have not been tested for MRI safety and compatibility.

OTHER INFORMATION

Documentation

For documentation purposes, remove the perforated section of the attached device label and attach it to the patient information card and implantation slip.

Information on materials and packaging

Device group	Devices	Material
Implants	BoneTrust® plus/sinus/cone+/cone+ sinus/cone PWR/one/mini/mini+ Implantate	Titanium Grade 4 (3.7065)
Abutments	Esthetic Abutments, ceramic base, wide body abutments, adhesive bases, double crown abutment, direct abutments, healing abutments	Titanium Grade 5 (3.7165)
	UCLA HSL abutments	Heraflat platinum/gold
	UCLA NEM abutments	Platinum/iridium
	"Lucky Lock" abutments, including 4-pce. Clip-Set	Abutment and male parts Titanium Grade 5 (3.7165) with TiN-coating Female parts: Nylon
Screws/caps	O-rings for BoneTrust mini	1.4305
	Healing caps, mounts and retaining screws	Titanium Grade 5 (3.7165)

Each pack consists of a glass container sealed tight with a plastic screw cap which contains the implant and the pre-attached mount. This container is sealed into a blister pack which acts as a sterilization barrier. The outer packaging consists of a labeled paper carton.

WARNING

The implant has been sterilized by radiation and must not be used after the expiration date. Nor should the implant be used if the packaging has been opened or is damaged.

CAUTION

The implants must be stored with their sealed packaging intact in a dry place at room temperature.

Expiration date/service life

Devices	Service life
BoneTrust® plus/sinus/cone+/cone+ sinus/cone PWR/ hex PWR/one	25 years
Abutments for BoneTrust® plus/sinus/cone+/cone+ sinus/cone PWR/hex PWR/one	25 years
BoneTrust® mini/mini+	9 months
Abutments for BoneTrust® mini/mini+	9 months

Accessories

Please refer to the company's homepage for the latest information on additional replacement parts and order numbers relating, for example to the placement tools.

www.medical-instinct.de

Disposal

Dispose of all devices and products used in accordance with the disposal guidelines for medical practices and hospitals.

MANUFACTURER

Medical Instinct® Deutschland GmbH
Graseweg 24
D-37120 Bovenden
www.medical-instinct.de

EXPLANATION OF SYMBOLS

Manufacturer

Expiration date

Follow the instruction or Use

Article number

Batch description

Radiation sterilized

Do not reuse

Do not re-sterilize

Do not use if the packaging is damaged

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– All rights reserved –
© Medical Instinct® Deutschland GmbH