

## Specifications



|                        |                       | MEDIT i700  | MEDIT i700 wireless   |
|------------------------|-----------------------|---|---|
| Scan-Technologie       | Scanbereich           | Up to 70 FPS  | Up to 70 FPS  |
|                        | Bildtechnologie       | 3D-in-motion video technology<br>3D full color streaming capture  | 3D-in-motion video technology<br>3D full color streaming capture  |
|                        | Leuchtmittel          | LED   | LED   |
|                        | Anti-Beschlag-Technik | Adaptive anti-fog   | Adaptive anti-fog   |
| Genauigkeit            | Komplettkiefer        | 10,9µm ± 0,98   | 10,9µm ± 0,98   |
| Handstück              | Maße                  | 248 x 44 x 47,4 mm  | 313 x 44 x 47,4 mm  |
|                        | Gewicht               | 245 g   | 328 g   |
| Scankopf               | Scankopfgröße         | 22,2 x 15,9 mm  | 22,2 x 17,1 mm  |
|                        | Spiegelwinkel         | 45 degree angle<br>Easier scanning of the distal molar region   | –   |
|                        | Scanfläche            | 15 x 13 mm  | 15 x 13 mm  |
|                        | Autoklavierbarkeit    | Bis zu 100 x<br>- autoclavable at 121 °C – 30 Min.<br>- autoclavable at 134 °C – 4 Min.<br>- autoclavable at 135 °C – 10 Min. | Bis zu 150 x<br>- autoclavable at 121 °C – 30 Min. (Gravity Type)<br>- autoclavable at 134 °C – 4 Min. (Pre-Vacuum Type)<br>- autoclavable at bei 135 °C – 10 Min. (Gravity Type) |
|                        | Drehbarer Scankopf    | Yes   | Yes   |
| Spezialmodus           | Fernsteuerungsmodus   | Yes   | Yes   |
|                        | UV-C Desinfektion     | Yes   | Yes   |
| Anschlussmöglichkeiten | Kabel / Wireless      | 2.0 m cable (optionally 2.5 m) / detachable handpiece   | Wireless  |
|                        | Anschluss             | USB 3.1 Gen1 (C Power Delivery)   | USB 3.1 Gen1 (Wireless Hub)   |