



## Order information

Scan-Caps and Ceramic Blocks (size "S") are available from your dealer / dental depot.



Scanbody S for Bluecam  
(36 pieces) - color white  
order number: 64 31 295



Scanbody S for Omnicam  
(36 pieces) - color grey  
order number: 64 31 311



Always use block size S.



## Order information

### BONETRUST® PLUS HEX IMPLANTS

■ ■ CAD/CAM Ti-Base hex for Cerec®,  
3.4/4.0 mm, with screw

■ CAD/CAM Ti-Base hex for Cerec®,  
5.0 mm, with screw

### SCANPOST HEX FOR CEREC®

■ ■ ■ Scanpost hex for CEREC®,  
3.4/4.0/5.0 mm, inkl. Schraube

### BONETRUST® PLUS CONE IMPLANTS

■ ■ CAD/CAM Ti-Base cone+ for Cerec®,  
3.4/4.0 mm, with screw

■ CAD/CAM Ti-Base cone+ for Cerec®,  
5.0 mm, with screw

### SCANPOST CONE FOR CEREC®

■ ■ Scanpost cone+ for CEREC®,  
3.4/4.0 mm, with screw

■ Scanpost cone+ for CEREC®,  
5.0 mm, inkl. Schraube

The BoneTrust® Titanium Bases are only available from Medical Instinct®

ART.-NO.:

166-303440

166-305000

166-323440

266-303440

266-305000

266-323440

266-325000

medical inst+inct®  
... WE ARE IMPLANTOLOGY!


Medical Instinct® Deutschland GmbH

Graseweg 24 . 37120 Bovenden, Germany  
Phone: 0 55 93 . 9 51 96 . E-Mail: info@medical-instinct.de

www.medical-instinct.de

# CAD/CAM TI - BASE for CEREC® DIGITAL WORKFLOW

Full integration into Sirona CEREC®\*  
with „Camlog®\*\* interface“



medical inst+inct®  
... WE ARE IMPLANTOLOGY!


\*Cerec is a registered trademark of Sirona Dental GmbH

\*\*Camlog® is a registered trademark of Camlog AG



# CEREC\* - Digital Workflow for Medical Instinct® Ti-Base

- Administration
- Define restoration
- Abutment




## ■ Scan plaster model or intra-oral

- Use Scanbody S for Bluecam (64 31 295)
  - Use Scanbody for Omnicam (64 31 311)
- (Scanbodies are available from your dealer / dental depot)



## ■ Material selection


- Sirona inCoris ZI meso
- Ivoclar Vivadent IPS e.max
- Always use blocks size "S"




## ■ Scanbodytype TiBase or Scanpost

### ■ Choose Camlog®

- Camlog® 3.8 mm for Medical Instinct® Ti-Base 3.4 / 4.0 mm
- Camlog® 4.3 mm for Medical Instinct® Ti-Base 5.0 mm



## ■ Abutment Design



## ■ Milling Abutment

